

# Zematra Total Sediment Tester

ASTM D4870; IP375; ISO 10307



**The Zematra 2 position Sediment Tester determines the insoluble material content of distillate and residual fuel oils.**

**It is a bench model which consists of two filtration cells, vacuum circuit with gauge and integrated cooling/heating circuit.**

#### **Optional:**

- Vacuum pump
- Steam generator
- Waste container
- Thermometer
- Balance

#### **Specifications:**

- Measures sediment up to 0.5% m/m
- Samples viscosities up to 55 mm<sup>2</sup>/s at 100°C
- 10 +/- 0.5 gram sample size
- 2 test stations
- Dimensions: 50x30x30cm.

**Each cup has to be equipped with a sintered brass disc to support the microfibre filters. This is a additional accessory.**



**Zematra, your reliable partner in testing**



2139 High Tech Road | State College, PA | 16803 | USA  
800-676-6232 | 814-353-8000 | Fax: 814-353-8007  
sales@cannoninstrument.com | cannoninstrument.com