



Automated Small Scale Closed Cup Flash Point Tester

asc-8c

Tanaka's 8th generation Flash Point Testers welcome a new member "asc-8c" for Small Scale Closed Cup Method.

Ergonomic

- Color coordinated display for high visibility and intuitive operation.
- Universal design

Robust

- Premium design and high quality finish for many years of reliable use.

Versatile

- Interchangeable ignition source from gas to electric or vice versa.
- Easily installable and removable maintenance parts for convenient use.



Specifications

Compliant Methods	ASTM D3278, D3828, D7236, ISO 3679, 3680, IP 523, 524 and others
Measuring Range (in ambient environment)	Default : -20 °C to 135 °C With external chiller : -30 °C to 135 °C
Ignition Source	Gas Flame or Electric Coil
Size (W x D x H)	230 x 470 x 390 (mm)
Weight	19 kg
Power Supply	AC100-240 V, 50/60 Hz, 400 VA

Specifications subject to change without prior notice.